beurer

BF 480 USB

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Dear Customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for the applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,
Your Beurer Team

1. Important notes – keep for later use!

1.1 Safety notes:

- **The scales must not be used by persons with medical implants (e.g. heart pace-makers). Otherwise their function could be impaired.**
- Do not use during pregnancy.
- Caution! Do not stand on the scale with wet feet and do not step on the scale when its surface is wet – danger of slipping!
- Batteries are highly dangerous if swallowed. Keep batteries and scale out of reach of small children. If batteries are swallowed, get medical help immediately.
- Keep children away from packing materials (danger of suffocation).
- Batteries must not be charged or reactivated by other means, taken apart, thrown into a fire or short-circuited.

1.2 General notes:

- Not intended for medical or commercial use.
- Please note that measuring tolerances are possible due to technical reasons because the scale has not been calibrated.
- Ages 10 ... 100 years and height settings from 100 ... 220 cm (3'-03"-7'-03") can be preset. Load capacity: max. 180 kg (396 lb, 28 st.). Result in steps of 100g (0.2 lb, 0.01st). Measurements of body fat, body water and muscle density in steps of 0.1%.
- When supplied to the customer, the scale is set to weigh and measure in "kg" and "cm". To the rear of the scales, there is a switch where you can select 'pound' and 'stone' (lb, st).
- Place the scale on a firm level floor; a firm floor covering is required for correct measurement.
- The scale should be cleaned occasionally with a damp cloth. Do not use abrasive detergents and never immerse the scale in water.
- Protect the scale against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stove, heating radiators).
- Repairs may only be carried out by customer service or authorized dealers. Before submitting any complaint, first check the batteries and replace them if necessary.

2. Operation

2.1 Inserting batteries

If present, pull the battery insulating strip off the battery compartment cover or remove the battery's protective film and insert the battery according to the polarity. If the scale fails to operate, remove the battery completely and insert it again. Your scale is equipped with a "replace battery" indicator. When the scale is operated with a battery which is too weak, "LU" appears in the display field and the scale automatically switches off. In this case the battery must be replaced (1 x 3V lithium battery CR2032).

2.2 Basic settings

Insert the batteries/"SET" (3 sec.)

Sequence	Set values
Date	year – month – day
Time	hours – minutes
Measurement unit	KG - LB - ST

Use the ◀ and ▶ buttons to select and confirm with "SET".

2.3 Setting user data

"SET"

User data	Set values
Memory position P	1 to 10
Height cm	100 to 220 cm (3'-03" to 7'-03")
Date of birth Age	year, month, day
Sex	male  , female 
Degree of activity A	1-5

Use the ◀ and ▶ buttons to select and confirm with "SET".

Degrees of activity

Degrees of activity	Physical activity
1	None.
2	Low: A small amount of light physical effort (e.g. short walks, light garden work, gymnastic exercises).
3	Medium: Physical effort for 30 minutes at least 2 to 4 times a week.
4	High: Physical effort for 30 minutes at least 4 to 6 times a week.
5	Very high: Intensive physical effort, intensive training or hard physical work for at least one hour daily.

3. Operation

3.1 Measuring weight only

- Step briefly and forcefully on the weighing area of your scale.
The entire display is shown (Fig. 1) until "0.0" (Fig. 2) is displayed.
- Stand on the scale without moving about and distribute your weight evenly on both legs. The scale immediately begins to measure your weight. Your measured weight is displayed shortly after this (Fig 3).
- The scale switches off a few seconds after you step down from it.

Fig. 1

Fig. 2

Fig. 3

3.2 Measuring weight, body fat, body water and muscle percentage

First your personal parameters must be stored;

- Step briefly on the weighing area of your scale.
- Select the memory preset where your basic personal data are stored by pressing the ► or ◀ button repeatedly. The data are displayed consecutively until "0.0" (Fig. 2) appears.
- Step onto the scale barefoot and make sure you're standing on both electrodes. First your body weight is determined and displayed (Fig. 3).
- Remain on the scale and avoid moving; your body fat is being analysed. This may take a few seconds.

-  When the unit is switched off, the measurements from an analysis are automatically saved. If the 60 memory locations for a user are full, "Full" appears on the display.

The following data are displayed:

- Body fat percentage BF in % (Fig. 4)
- Water percentage ≈ in % (Fig. 5)
- Muscle percentage ↔ in % (Fig. 6)
- Now all measured values are displayed consecutively and the scale switches off.
- Weight in kg (lb, st)

Fig. 4


Fig. 5

Fig. 6

Fig. 3

3.3 After a body fat measurement

After the measurement is complete, the display reads "0.0".

- Press the ◀ button. The memory symbol M appears with the most recently stored measurement. The results for body fat, body water, muscle percentage and weight are displayed in succession.
- Press the ◀ button again. The second most recently saved measurement appears.
- Press the ◀ button again. The memory symbol goes out and "0.0" appears in the display.

-  To delete the measured value shown, press the "SET" button, and confirm and delete "dE!S" by pressing the "SET" button again.
- Clearing the memory: After switching on, press and hold the ◀, ► and "SET" buttons simultaneously for 5 seconds until "CLL" appears. Confirm "YES" and delete by pressing the "SET" button.

3.4 Tips for using the scale

Important when measuring body fat/body water/muscle percentage:

- The measurement may only be carried out while barefoot and it is helpful if the soles of your feet are slightly damp. Completely dry soles can result in unsatisfactory results, as they have insufficient conductivity.
- Your feet, legs, calves and thighs must not touch each other, as otherwise the measurement cannot be performed correctly.
- Stand still during the measurement.
- Wait several hours (6-8) after unusually strenuous activity.
- Wait approx. 15 minutes after getting out of bed so that the water in your body can be distributed.

The measurement is not reliable for:

- Children under approx. 10 years of age.
- Persons with fever, undergoing dialysis, with symptoms of edema or osteoporosis.
- Persons taking cardiovascular medication.

- Persons taking vasodilating or vasoconstricting medications.
- Persons with substantial anatomical deviations in the legs relative to their total height (leg length considerably shorter or longer than usual).

4. Evaluation of results

Body fat percentage

The following body fat percentages are guide values (contact your physician for further information).

Man

Age	low	normal	high	very high
10-14	<11	11-16	16,1-21	>21,1
15-19	<12	12-17	17,1-22	>22,1
20-29	<13	13-18	18,1-23	>23,1
30-39	<14	14-19	19,1-24	>24,1
40-49	<15	15-20	20,1-25	>25,1
50-59	<16	16-21	21,1-26	>26,1
60-69	<17	17-22	22,1-27	>27,1
70-100	<18	18-23	23,1-28	>28,1

Woman

Age	low	normal	high	very high
10-14	<16	16-21	21,1-26	>26,1
15-19	<17	17-22	22,1-27	>27,1
20-29	<18	18-23	23,1-28	>28,1
30-39	<19	19-24	24,1-29	>29,1
40-49	<20	20-25	25,1-30	>30,1
50-59	<21	21-26	26,1-31	>31,1
60-69	<22	22-27	27,1-32	>32,1
70-100	<23	23-28	28,1-33	>33,1

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values stated. It should, however, be noted that there could be a danger to health in the case of extremely low values.

Body water percentage

The body water percentage is normally within the following ranges:

Man

Age	poor	good	very good
10-100	<50	50-65	>65

Woman

Age	low	normal	high
10-100	<65	45-60	>60

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age-related water retention, for example. If necessary ask your physician. Basically, a high body water percentage should be the aim.

Muscle percentage

The muscle percentage % is normally within the following ranges:

Man

Age	low	normal	high
10-14	<44	44-57	>57
15-19	<43	43-56	>56
20-29	<42	42-54	>54
30-39	<41	41-52	>52
40-49	<40	40-50	>50
50-59	<39	39-48	>48
60-69	<38	38-47	>47
70-100	<37	37-46	>46

Woman

Age	low	normal	high
10-14	<36	36-43	>43
15-19	<35	35-41	>41
20-29	<34	34-39	>39
30-39	<33	33-38	>38
40-49	<31	31-36	>36
50-59	<29	29-34	>34
60-69	<28	28-33	>33
70-100	<27	27-32	>32

5. PC transfer


System requirements: Windows 7, Windows Vista, Windows XP SP3

Installation process:

- You will find the Beurer software download and installation in the download area at <http://www.beurer.com>
- **HealthManager:** Health management in conjunction with Beurer blood pressure and blood glucose monitors.
- Follow the instructions for software installation.

Transfer data:

- Establish USB connection.
- Select data transfer in the software.

6. Incorrect measurement

If the scale detects an error during weighing, "FFFF" or "Err." appears in the display.

If you step onto the scale before "0.0" appears in the display, the scale will not operate properly.

Possible causes of errors:	Remedy:
– The maximum load-bearing capacity of 180 kg (396 lbs) was exceeded.	– Only weigh the maximum permissible weight.
– The electrical resistance between the electrodes and the soles of your feet is too high (e.g. with heavily callused skin).	– Repeat weighing barefoot. – Slightly moisten the soles of your feet if necessary. Remove the calluses on the soles of your feet if necessary.
– Your body fat lies outside the measurable range (less than 5% or greater than 55%).	– Repeat weighing barefoot. – Slightly moisten the soles of your feet if necessary.
– Scale is on carpeting.	– Place the scale on a level and solid surface.

7. Disposal

Standard and rechargeable batteries should not be disposed of separately from the household waste. As a consumer, you are legally obliged to return used batteries for proper disposal. You can hand in your used batteries at public collection points in your district or sales outlets where batteries of this type are sold.

Note

The codes below are printed on batteries containing harmful substances: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

For environmental reasons, do not dispose of the scale with the batteries in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Observe the local regulations for material disposal. Dispose of the unit in accordance with EC Directive 2002/96/EC– WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

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